

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Adenosine Kinase (ADK) ELISA KBS-KBH6011

|                            |                                                                                                                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Adenosine Kinase (ADK) ELISA                                                                                                                                                                                                                          |
| Biozol Catalog Number      | KBS-KBH6011                                                                                                                                                                                                                                                         |
| Supplier Catalog Number    | KBH6011                                                                                                                                                                                                                                                             |
| Alternative Catalog Number | KBS-KBH6011-1                                                                                                                                                                                                                                                       |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                                 |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                         |
| Species Reactivity         | Human                                                                                                                                                                                                                                                               |
| Product Description        | The GENLISA Human Adenosine Kinase (ADK) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use<br>Standardisation and High Reproducibility Lot to Lot Consistency<br>Accuracy and Precision Validated against seven points ... |
| Range                      | 1.56 - 100 ng/ml                                                                                                                                                                                                                                                    |
| Sensitivity                | 0.75 ng/ml                                                                                                                                                                                                                                                          |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                                  |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                          |