

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Cathepsin L3 (CTSL3) ELISA KBS-KBH9663

|                            |                                                                                                                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | GENLISA Human Cathepsin L3 (CTSL3) ELISA                                                                                                                                                                                                                            |
| Biozol Catalog Number      | KBS-KBH9663                                                                                                                                                                                                                                                         |
| Supplier Catalog Number    | KBH9663                                                                                                                                                                                                                                                             |
| Alternative Catalog Number | KBS-KBH9663-1                                                                                                                                                                                                                                                       |
| Manufacturer               | Krishgen Biosystems                                                                                                                                                                                                                                                 |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                         |
| Species Reactivity         | Human                                                                                                                                                                                                                                                               |
| Product Description        | The GENLISA Human Cathepsin L3 (CTSL3) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use<br>Standardisation and High Reproducibility Lot to Lot Consistency<br>Accuracy and Precision Validated against seven points fo... |
| Range                      | 0 - 1000 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                          |
| Sensitivity                | 30 ng/ml* (May vary. Please confirm)                                                                                                                                                                                                                                |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                                  |
| Application Notes          | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                          |