

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## **GENLISA Human Filaggrin 2 (FLG2) ELISA KBS-KBH9855**

Article Name	GENLISA Human Filaggrin 2 (FLG2) ELISA
Biozol Catalog Number	KBS-KBH9855
Supplier Catalog Number	KBH9855
Alternative Catalog Number	KBS-KBH9855-1
Manufacturer	Krishgen Biosystems
Category	Kits/Assays
Species Reactivity	Human
Product Description	The GENLISA Human Filaggrin 2 (FLG2) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precision Validated against seven points for
Range	0 - 1000 ng/ml* (May vary. Please confirm)
Sensitivity	30 ng/ml* (May vary. Please confirm)
Samples	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.
Application Notes	Colorimetric, Absorbance measured at 450nm