

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### GENLISA Human Histone Deacetylase 10 (HDAC10) ELISA KBS-KBH9945

|                            |   |
|----------------------------|---|
| Article Name               | GENLISA Human Histone Deacetylase 10 (HDAC10) ELISA   |
| Biozol Catalog Number      | KBS-KBH9945   |
| Supplier Catalog Number    | KBH9945   |
| Alternative Catalog Number | KBS-KBH9945-1   |
| Manufacturer               | Krishgen Biosystems   |
| Category                   | Kits/Assays   |
| Species Reactivity         | Human   |
| Product Description        | The GENLISA Human Histone Deacetylase 10 (HDAC10) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precision Validated against seve... |
| Range                      | 0 - 1000 ng/ml* (May vary. Please confirm)  |
| Sensitivity                | 30 ng/ml* (May vary. Please confirm)  |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.  |
| Application Notes          | Colorimetric, Absorbance measured at 450nm  |