

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## Recombinant Human Small ubiquitin-related modifier 5 (SUMO5) Preis auf Anfrage CSB-BP140323HU

Article Name	Recombinant Human Small ubiquitin-related modifier 5 (SUMO5) Preis auf Anfrage
Biozol Catalog Number	CSB-BP140323HU
Supplier Catalog Number	CSB-BP140323HU
Alternative Catalog Number	CSB-BP140323HU-20, CSB-BP140323HU-100, CSB-BP140323HU-1
Manufacturer	Cusabio
Category	Proteine/Peptide
Species Reactivity	Human
Molecular Weight	1
Tag	Tag type will be determined during the manufacturing process.  (Note:If you have specified tag type, please tell us or remark on your PO and we will develop the specified tag preferentially.)
UniProt	G2XKQ0
Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 20%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Purity	Greater than 85% as determined by SDS-PAGE.
Form	Liquid or Lyophilized powder
Sequence	Full Length of Mature Protein

Target	MSDLEAKPSTEHLGDKIKDEDIKLRVIGQDSSEIHFKVKMTTPLKKLKKSYCQR QGVPVNSLRFLFEGQRIADNHTPEELGMEEEDVIEVYQEQIGG
Application Notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.