

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Adhesion Molecule Array Q1 RAY-QAH-CAM-1-2

|                            |                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human Adhesion Molecule Array Q1                                                                                                   |
| Biozol Catalog Number      | RAY-QAH-CAM-1-2                                                                                                                    |
| Supplier Catalog Number    | QAH-CAM-1-2                                                                                                                        |
| Alternative Catalog Number | RAY-QAH-CAM-1-2                                                                                                                    |
| Manufacturer               | RayBiotech                                                                                                                         |
| Category                   | Kits/Assays                                                                                                                        |
| Species Reactivity         | Human                                                                                                                              |
| Product Description        | Quantibody Human Adhesion Molecule Array 1 Kit. Detects 17 human cell adhesion molecules. Suitable for all liquid sample types.... |
| Samples                    | Cell Culture Supernatants, Plasma, Serum, Tissue Lysates, Other Body Fluids, Cell Lysates                                          |
| Target                     | VE-Cadherin (CDH5)                                                                                                                 |