

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human Ig Isotyping Array Q1 RAY-QAH-ISO-1-1

|                            |                                                                                                                                 |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human Ig Isotyping Array Q1                                                                                                     |
| Biozol Catalog Number      | RAY-QAH-ISO-1-1                                                                                                                 |
| Supplier Catalog Number    | QAH-ISO-1-1                                                                                                                     |
| Alternative Catalog Number | RAY-QAH-ISO-1-1                                                                                                                 |
| Manufacturer               | RayBiotech                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                     |
| Species Reactivity         | Human                                                                                                                           |
| Product Description        | Quantibody Human Immunoglobulin Isotyping Array Kit. Detects 8 Human Immunoglobulin Sub-Types. Suitable for serum or plasma.... |
| Samples                    | Cell Culture Supernatants, Plasma, Serum, Tissue Lysates, Other Body Fluids, Cell Lysates                                       |
| Target                     | IgA, IgD, IgE, IgG1, IgG2, IgG3, IgG4, IgM                                                                                      |