

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Custom Mouse TIMP-1 SpeedELISA Kit RAY-ELMS-TIMP1-2

|                            |                                                                                                                                                                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Custom Mouse TIMP-1 SpeedELISA Kit                                                                                                                                                                                                                      |
| Biozol Catalog Number      | RAY-ELMS-TIMP1-2                                                                                                                                                                                                                                        |
| Supplier Catalog Number    | ELMS-TIMP1-2                                                                                                                                                                                                                                            |
| Alternative Catalog Number | RAY-ELMS-TIMP1-2                                                                                                                                                                                                                                        |
| Manufacturer               | RayBiotech                                                                                                                                                                                                                                              |
| Category                   | Kits/Assays                                                                                                                                                                                                                                             |
| Species Reactivity         | Mouse                                                                                                                                                                                                                                                   |
| Product Description        | RayBio Mouse Custom TIMP-1 SpeedELISA is a sandwich-based immunoassay with a simplified work flow for the rapid, quantitative measurement of mouse TIMP-1 in only 3 hours. Suitable for serum, plasma, cell culture supernatants, or lysate samples.... |
| NCBI                       | <a href="#">21857</a>                                                                                                                                                                                                                                   |
| Form                       | 96-well Microplate                                                                                                                                                                                                                                      |
| Target                     | TIMP-1 (TIMP1)                                                                                                                                                                                                                                          |
| Application Notes          | Colorimetric                                                                                                                                                                                                                                            |