

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

S-Nitrosylation Detection Kit (Biotin Switch) RAY-MA-SNO-W

| | |
|----------------------------|--|
| Article Name | S-Nitrosylation Detection Kit (Biotin Switch) |
| Biozol Catalog Number | RAY-MA-SNO-W |
| Supplier Catalog Number | MA-SNO-W |
| Alternative Catalog Number | RAY-MA-SNO-W |
| Manufacturer | RayBiotech |
| Category | Kits/Assays |
| Product Description | RayBio S-Nitrosylation Detection Kit uses a modified biotin-switch method to allow for the direct visualization of S-nitrosylated proteins by western blot analysis. This assay can label 40 samples containing 100-200 µg of total protein each.... |