

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### DNA Bisulfite Conversion Kit (Magnetic Bead Method) RAY-NAE-BIS-MAG

Article Name	DNA Bisulfite Conversion Kit (Magnetic Bead Method)
Biozol Catalog Number	RAY-NAE-BIS-MAG
Supplier Catalog Number	NAE-BIS-MAG
Alternative Catalog Number	RAY-NAE-BIS-MAG
Manufacturer	RayBiotech
Category	Kits/Assays
Species Reactivity	Human
Product Description	This kit is used to convert unmethylated cytosine residues into uracil for gene methylation analysis in genomic DNA. Methylated cytosine residues (5-methylcytosine) remain unchanged. PCR primers and/or probes can be designed to distinguish methylated...
Samples	DNA