

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

Complement 4d (Acidic Complement C4, Basic C4, C4, C4-1, C4a Anaphylatoxin, C4A, C4A2, C4A3, C4A4, C4A6, C4b, C4d, C4B1, C4B12, C4d fragment, C4F, C4S, Complement C4 Gamma Chain, Complement Component 4A, Complement Component 4B) ( USB-172070-AF

| Article Name               | Complement 4d (Acidic Complement C4, Basic C4, C4, C4-1, C4a<br>Anaphylatoxin, C4A, C4A2, C4A3, C4A4, C4A6, C4b, C4d, C4B1,<br>C4B12, C4d fragment, C4F, C4S, Complement C4 Gamma Chain,<br>Complement Component 4A, Complement Component 4B) (            |
|----------------------------|--|
| Biozol Catalog Number      | USB-172070-AF  |
| Supplier Catalog Number    | 172070-AF  |
| Alternative Catalog Number | USB-172070-AF-100  |
| Manufacturer               | US Biological  |
| Host                       | Mouse  |
| Category                   | Antikörper   |
| Application                | IF, IHC  |
| Immunogen                  | Recombinant human Complement 4d protein. Cellular Localization:<br>Intracytoplasmic vacuoles of endothelial cells and Secreted.  |
| Product Description        | This MAb is specific to Complement 4d (C4d) and it reacts with the secreted as well as cell-bound C4d. C4d is a degradation product of the activated complement factor C4b. Complement 4b is typically activated by binding of Abs to specific target mole |
| Clonality                  | Monoclonal   |
| Clone Designation          | [SPM545]   |

| UniProt | P0C0L4  |
|---------|---|
| Purity  | Purified by Protein A/G affinity chromatography.                            |
| Form    | Supplied as a liquid in PBS. Also available with BSA and azide, see 172070. |