

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Breast Cancer Resistance Protein (BCRP, BCRP1, ABC15, ABCP, ATP-binding Cassette Sub-family G Member 2, ABCG2, BMDP, CD338, CDw338, EST157481, MGC102821, Mitoxantrone Resistance-associated Protein, MRX, MXR, MXR1, Placenta-specific 1, P-glycoprotein, P-gp, P-gp1, P-gp2, P-gp3, P-gp4, P-gp5, P-gp6, P-gp7, P-gp8, P-gp9, P-gp10, P-gp11, P-gp12, P-gp13, P-gp14, P-gp15, P-gp16, P-gp17, P-gp18, P-gp19, P-gp20, P-gp21, P-gp22, P-gp23, P-gp24, P-gp25, P-gp26, P-gp27, P-gp28, P-gp29, P-gp30, P-gp31, P-gp32, P-gp33, P-gp34, P-gp35, P-gp36, P-gp37, P-gp38, P-gp39, P-gp40, P-gp41, P-gp42, P-gp43, P-gp44, P-gp45, P-gp46, P-gp47, P-gp48, P-gp49, P-gp50, P-gp51, P-gp52, P-gp53, P-gp54, P-gp55, P-gp56, P-gp57, P-gp58, P-gp59, P-gp60, P-gp61, P-gp62, P-gp63, P-gp64, P-gp65, P-gp66, P-gp67, P-gp68, P-gp69, P-gp70, P-gp71, P-gp72, P-gp73, P-gp74, P-gp75, P-gp76, P-gp77, P-gp78, P-gp79, P-gp80, P-gp81, P-gp82, P-gp83, P-gp84, P-gp85, P-gp86, P-gp87, P-gp88, P-gp89, P-gp90, P-gp91, P-gp92, P-gp93, P-gp94, P-gp95, P-gp96, P-gp97, P-gp98, P-gp99, P-gp100)** **USB-B2738-01D**

Article Name	Breast Cancer Resistance Protein (BCRP, BCRP1, ABC15, ABCP, ATP-binding Cassette Sub-family G Member 2, ABCG2, BMDP, CD338, CDw338, EST157481, MGC102821, Mitoxantrone Resistance-associated Protein, MRX, MXR, MXR1, Placenta-specific 1, P-glycoprotein, P-gp, P-gp1, P-gp2, P-gp3, P-gp4, P-gp5, P-gp6, P-gp7, P-gp8, P-gp9, P-gp10, P-gp11, P-gp12, P-gp13, P-gp14, P-gp15, P-gp16, P-gp17, P-gp18, P-gp19, P-gp20, P-gp21, P-gp22, P-gp23, P-gp24, P-gp25, P-gp26, P-gp27, P-gp28, P-gp29, P-gp30, P-gp31, P-gp32, P-gp33, P-gp34, P-gp35, P-gp36, P-gp37, P-gp38, P-gp39, P-gp40, P-gp41, P-gp42, P-gp43, P-gp44, P-gp45, P-gp46, P-gp47, P-gp48, P-gp49, P-gp50, P-gp51, P-gp52, P-gp53, P-gp54, P-gp55, P-gp56, P-gp57, P-gp58, P-gp59, P-gp60, P-gp61, P-gp62, P-gp63, P-gp64, P-gp65, P-gp66, P-gp67, P-gp68, P-gp69, P-gp70, P-gp71, P-gp72, P-gp73, P-gp74, P-gp75, P-gp76, P-gp77, P-gp78, P-gp79, P-gp80, P-gp81, P-gp82, P-gp83, P-gp84, P-gp85, P-gp86, P-gp87, P-gp88, P-gp89, P-gp90, P-gp91, P-gp92, P-gp93, P-gp94, P-gp95, P-gp96, P-gp97, P-gp98, P-gp99, P-gp100)
Biozol Catalog Number	USB-B2738-01D
Supplier Catalog Number	B2738-01D
Alternative Catalog Number	USB-B2738-01D-1
Manufacturer	US Biological
Host	Mouse
Category	Antikörper
Application	FC, IHC, WB
Immunogen	Mitoxanthrone resistant, BCRP overexpressing cell line MCF7 MR.
Product Description	Applications: Suitable for use in Flow Cytometry, Immunohistochemistry and Immunocytochemistry. Not suitable for use in Western Blot. Other applications not tested. Recommended Dilution: Flow Cytometry: 1:20. Permeabilisation necessary. Immunohistochem...
Clonality	Monoclonal
Clone Designation	[BXP-34]

Isotype	IgG1
Purity	Ascites
Form	Supplied as a liquid, 1% BSA, 0.09% sodium azide.