



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Glycine max (SBA) (TRITC) USB-G8166-30

Article Name	Glycine max (SBA) (TRITC)
Biozol Catalog Number	USB-G8166-30
Supplier Catalog Number	G8166-30
Alternative Catalog Number	USB-G8166-30-2
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Purified Glycine max (SBA) covalently linked to agarose beads (exclusion limit MW 5-15 million). Lectins are derived from the extracts of plants, animals, viruses and microorganisms and are known to agglutinate red blood cells. These agglutinins can
Form	Supplied as a liquid in 0.01M phosphate, 0.15M sodium chloride, pH 7.2-7.4, 0.05% sodium azide. Labeled with TRITC.