



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

IgG Fab, Equine (FITC) USB-I1903-20Q

Article Name	IgG Fab, Equine (FITC)
Biozol Catalog Number	USB-I1903-20Q
Supplier Catalog Number	I1903-20Q
Alternative Catalog Number	USB-I1903-20Q-1
Manufacturer	US Biological
Category	Molekularbiologie
Application	FLISA, IF, WB
Product Description	Label: Fluorescein isothiocyanate (FITC) (Molecular Weight 390 daltons) Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 1.6 moles FITC per mole of Horse IgG Fab Applications: Suitable for use in FLISA (Fluorescen
Purity	Purified by ion exchange chromatography. IEP Assay shows a single precipitin arc against source serum protein and IgG.
Form	Supplied as a lyophilized powder in PBS, pH 7.2, 10mg/ml BSA, 0.09% sodium azide. Labeled with Fluorescein isothiocyanate (FITC).