



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

IgG, Feline (FITC) USB-I1903-22N

Article Name	IgG, Feline (FITC)
Biozol Catalog Number	USB-I1903-22N
Supplier Catalog Number	I1903-22N
Alternative Catalog Number	USB-I1903-22N-1
Manufacturer	US Biological
Category	Molekularbiologie
Application	FLISA, IF
Product Description	Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Fluorochrome/Protein Ratio: 2.1 moles FITC per mole of Cat IgG Applications: Suitable for use in FLISA and Immunofluorescence. Other applications have not been tested. Recommended Dilutions: O
Purity	Purified by ion exchange chromatography. IEP Assay shows a single precipitin arc against source serum protein and IgG.
Form	Supplied as a lyophilized powder in PBS, pH 7.2, 10mg/ml BSA, 0.09% sodium azide. Labeled with Fluorescein isothiocyanate (FITC). Reconstitute with 1ml sterile ddH2O.