



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

**Telephone:** +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## **Product Datasheet**

## IgG Fab, Goat (FITC) USB-I1903-25M

Article Name	IgG Fab, Goat (FITC)
Biozol Catalog Number	USB-I1903-25M
Supplier Catalog Number	I1903-25M
Alternative Catalog Number	USB-I1903-25M-1
Manufacturer	US Biological
Category	Molekularbiologie
Application	FLISA, IF, WB
Product Description	Absorption Wavelength: 495 nm Emission Wavelength: 528 nm Applications: Suitable for use in FLISA and Immunofluorescence as well as other antibody based fluorescent assays. Other applications have not been tested. Storage and Stability: Lyophilized p
Purity	Purified by ion exchange chromatography. IEP Assay shows a single precipitin arc against source serum protein and IgG.
Form	Supplied as a lyophilized powder in PBS, pH 7.2, 10mg/ml BSA (Immunoglobulin and Protease free), 0.01% sodium azide. Reconstitute with 1ml sterile ddH2O. Labeled with Fluorescein isothiocyanate (FITC).