



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

IgG, Goat (Agarose) USB-I1903-25R

Article Name	IgG, Goat (Agarose)
Biozol Catalog Number	USB-I1903-25R
Supplier Catalog Number	I1903-25R
Alternative Catalog Number	USB-I1903-25R-10
Manufacturer	US Biological
Category	Molekularbiologie
Product Description	Applications: Suitable for use in antibody purification or to remove unwanted crossreactivity. Other applications not tested. Recommended Dilution: Optimal dilutions to be determined by the researcher. Storage and Stability: Stored at 4C. Do not free
Purity	Purified by ion exchange chromatography. IEP Assay shows a single precipitin arc against anti-goat serum prior to coupling to the beads. IgG fraction of antibody is coupled to activated agarose.
Form	Supplied as a liquid in PBS, pH 7.2, 0.05% sodium azide. No stabilizing proteins added. IgG fraction of antibody is coupled to activated agarose.