



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mueller Hinton Medium (Agar) (Powder) USB-M5000

Article Name	Mueller Hinton Medium (Agar) (Powder)
Biozol Catalog Number	USB-M5000
Supplier Catalog Number	M5000
Alternative Catalog Number	USB-M5000-500
Manufacturer	US Biological
Category	Biochemikalien
Product Description	Mueller Hinton Agar is used in antimicrobial susceptibility testing by the disk diffusion method. This formula conforms to National Committee for Clinical Laboratory Standards (NCCLS). Quality Control: Tested according to USB Standard Test Methods an