



Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | **Fax:** +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

TMB X-tra-Prestained KEM-4810H

Article Name	TMB X-tra-Prestained
Biozol Catalog Number	KEM-4810H
Supplier Catalog Number	4810H
Alternative Catalog Number	KEM-4810H-100ML
Manufacturer	Kementec
Category	Sonstiges
Application	ELISA
Form	Ready-to-use (RTU)
Storage	2-8°C in the dark
Application Notes	HRP in ELISA – high sensitivity - stained pink