

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

ActiveMax DLL4 muBeads, premium grade (for cells, 20 mum) ABS-MBS-C024-2.5MG

Article Name	ActiveMax DLL4 muBeads, premium grade (for cells, 20 mum)
Biozol Catalog Number	ABS-MBS-C024-2.5MG
Supplier Catalog Number	MBS-C024-2.5mg
Alternative Catalog Number	ABS-MBS-C024-2.5MG
Manufacturer	AcroBiosystems
Category	Sonstiges
Product Description	ActiveMax DLL4 μ Beads, premium grade (for cells, 20 μ m) - Sonstiges
Form	Powder
Target	DLL4
Application Notes	ActiveMax Human DLL4 μ Beads, premium grade (for cells, 20 μ m) is designed to synergistically activate Notch signaling in hematopoietic stem/progenitor cells (HSPCs). The synergistic interactions between Notch ligand Delta-like 4 and VCAM-1 are leveraged to enhance Notch signaling to support efficient T-cell lineage commitment, so the beads is recommended to use for in vitro production of induced pluripotent stem cell derived T cells.