

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### TNFSF18 Antibody (Biotin), Rabbit, Polyclonal BYT-ORB1272726

|                            |  |
|----------------------------|--|
| Article Name               | TNFSF18 Antibody (Biotin), Rabbit, Polyclonal  |
| Biozol Catalog Number      | BYT-ORB1272726   |
| Supplier Catalog Number    | orb1272726   |
| Alternative Catalog Number | BYT-ORB1272726-0.05  |
| Manufacturer               | Biorbyt  |
| Host                       | Rabbit   |
| Category                   | Antikörper   |
| Application                | ELISA, WB  |
| Species Reactivity         | Human  |
| Immunogen                  | Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hAITRL. Human AITRL specific antibody was purified by affinity chromatography and then biotinylated. |
| Conjugation                | Biotin   |
| Product Description        | TNFSF18 Antibody (Biotin)...   |
| Clonality                  | Polyclonal   |
| Concentration              | batch dependent  |
| NCBI                       | <a href="#">9</a>  |
| UniProt                    | <a href="#">Q9UNG2</a>   |
| Form                       | Lyophilized  |

|                   |  |
|-------------------|--|
| Target            | TNFSF18  |
| Application Notes | <p>Application Notes: ELISA:Sandwich:To detect hAITRL by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.25 - 1.0 µg/mL of this antibody is required. This biotinylated polyclonal antibody, in conjunction with our polyclonal Anti-Human AITRL as a capture antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hAITRL. Western Blot:To detect hAITRL by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant hAITRL is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions</p> |