

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Human STEAP2 full length protein-synthetic nanodisc BYT-ORB1743228

|                            |                                                                                                                                                                                                                                                                 |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Human STEAP2 full length protein-synthetic nanodisc                                                                                                                                                                                                             |
| Biozol Catalog Number      | BYT-ORB1743228                                                                                                                                                                                                                                                  |
| Supplier Catalog Number    | orb1743228                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | BYT-ORB1743228-10,BYT-ORB1743228-100,BYT-ORB1743228-50                                                                                                                                                                                                          |
| Manufacturer               | Biorbyt                                                                                                                                                                                                                                                         |
| Category                   | Proteine/Peptide                                                                                                                                                                                                                                                |
| Species Reactivity         | Human                                                                                                                                                                                                                                                           |
| Product Description        | Human STEAP2 full length protein-synthetic nanodisc...                                                                                                                                                                                                          |
| Molecular Weight           | The human full length STEAP2 protein has a MW of 56.1 kDa                                                                                                                                                                                                       |
| UniProt                    | <a href="#">Q8NFT2</a>                                                                                                                                                                                                                                          |
| Form                       | Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with pH lowe |
| Target                     | STEAP2                                                                                                                                                                                                                                                          |