

Diagnostica Vertrieb GmbH, Oehleckerring 11-13

22419 Hamburg, Germany

Telephone: +49 (0)89 3799666-6 | Fax: +49 (0)89 3799666-99

E-Mail: info@biozol.de

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse Matrix Metalloproteinase 7 (MMP-7) ELISA Kit BYT-ORB1807857

Article Name	Mouse Matrix Metalloproteinase 7 (MMP-7) ELISA Kit
Biozol Catalog Number	BYT-ORB1807857
Supplier Catalog Number	orb1807857
Alternative Catalog Number	BYT-ORB1807857-48, BYT-ORB1807857-96
Manufacturer	Biorbyt
Category	Kits/Assays
Species Reactivity	Mouse
Product Description	Mouse Matrix Metalloproteinase 7 (MMP-7) ELISA Kit
Range	62.5-4000pg/mL
Sensitivity	37.50 pg/mL
UniProt	Q10738
Samples	serum, plasma, Tissue homogenate and Other biological samples
Target	MMP-7

Application Notes

This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse MMP-7. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Mouse MMP-7 and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse MMP-7, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm 2 nm. The OD value is proportional to the concentration of Mouse MMP-7. You can calculate the concentration of Mouse MMP-7 in the samples by comparing the OD of the samples to the standard curve.