

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Bovine GH (Growth Hormone) ELISA Kit EKL-EBV27430

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Bovine GH (Growth Hormone) ELISA Kit                                                                                                                                                                                                                          |
| Biozol Catalog Number      | EKL-EBV27430                                                                                                                                                                                                                                                  |
| Supplier Catalog Number    | EBv27430                                                                                                                                                                                                                                                      |
| Alternative Catalog Number | EKL-EBV27430-96T                                                                                                                                                                                                                                              |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Bovine                                                                                                                                                                                                                                                        |
| Product Description        | This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Bovine GH. During the reaction, Bovine GH in the sample or standard competes with a fixed amount of Bovine GH on the solid phase ... |
| Range                      | 123.5-10000pg/mL                                                                                                                                                                                                                                              |
| Sensitivity                | 50.2pg/mL                                                                                                                                                                                                                                                     |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | GH                                                                                                                                                                                                                                                            |