

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rat EAAT2 (Excitatory Amino Acid Transporter 2) ELISA Kit EKL-ER21503

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | Rat EAAT2 (Excitatory Amino Acid Transporter 2) ELISA Kit                                                                                                                                                                                                     |
| Biozol Catalog Number      | EKL-ER21503                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | ER21503                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | EKL-ER21503-96T                                                                                                                                                                                                                                               |
| Manufacturer               | EnkiLife                                                                                                                                                                                                                                                      |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Species Reactivity         | Rat                                                                                                                                                                                                                                                           |
| Product Description        | This ELISA kit uses the Sandwich-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Rat EAAT2. Standards or samples are added to the micro ELISA plate wells and combined with the specific anti... |
| Range                      | 0.15-10ng/mL                                                                                                                                                                                                                                                  |
| Sensitivity                | 0.054ng/mL                                                                                                                                                                                                                                                    |
| Samples                    | Serum,plasma and other biological fluids                                                                                                                                                                                                                      |
| Target                     | EAAT2                                                                                                                                                                                                                                                         |