

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **SARS-CoV-2 Neutralizing Antibody ELISA Kit ABX392372**

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| Article Name               | SARS-CoV-2 Neutralizing Antibody ELISA Kit                                                              |
| Biozol Catalog Number      | ABX392372                                                                                               |
| Supplier Catalog Number    | abx392372                                                                                               |
| Alternative Catalog Number | ABX392372-96TESTS                                                                                       |
| Manufacturer               | Abbexa                                                                                                  |
| Category                   | Kits/Assays                                                                                             |
| Application                | ELISA                                                                                                   |
| Species Reactivity         | All                                                                                                     |
| Range                      | 1.563 ng/ml - 100 ng/ml                                                                                 |
| Sensitivity                | < 0.938 ng/ml                                                                                           |
| Samples                    | Serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Target                     | SARS-CoV-2 Neutralizing Antibody                                                                        |