

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Protein Molecular Weight Marker, Unconjugated Preis auf Anfrage BYT-ORB1784560

Article Name	Protein Molecular Weight Marker, Unconjugated Preis auf Anfrage
Biozol Catalog Number	BYT-ORB1784560
Supplier Catalog Number	orb1784560
Alternative Catalog Number	BYT-ORB1784560-5
Manufacturer	Biorbyt
Category	Proteine/Peptide
Application	SDS-PAGE
Conjugation	Unconjugated
Product Description	Protein Molecular Weight Marker...
Concentration	1X
Buffer	Preservative: 0.01% (w/v) Sodium Azide. Stabilizer: None. 2.5 mM Tris-HCl, pH 6.8, 2% (w/v) SDS, 250 mM NaCl, 1mM EDTA, 4% (v/v) beta-Mercaptoethanol, 35% (v/v) Glycerol and 0.05% (w/v) bromophenol blue
Purity	Protein Molecular Weight Marker is supplied in 1X loading buffer and is ready-to-use. It contains: [62.5 mM Tris-HCl, pH 6.8, 2% (w/v) SDS, 250 mM NaCl, 1mM EDTA, 4% (v/v) beta-Mercaptoethanol, 35% (v/v) Glycerol and 0.05% (w/v) bromophenol blue.] This product are unstained markers optimized for use with Coomassie blue staining, but the loading volume can be adjusted for compatibility with other staining methods (i.e. silver staining). The intensity of the 67 kDa bovine albumin band has been increased to serve as a reference.
Form	Liquid (in 1x Loading Buffer)

Application Notes	Application Notes: Protein Molecular Weight Marker has been tested and is designed for use in accurate estimation of protein molecular weight in SDS-polyacrylamide gel electrophoresis and western blotting. Thaw the cocktail at room temperature. Gently invert vial several times to ensure that the solution is homogeneous and that any precipitated material is re-dissolved. A loading volume of 5-10μl is appropriate for most gel formats
-------------------	---