

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Protein A-coupled Magnetic Beads (recommended for MPCLIA), E. coli ABS-MPC-A001-10MG

Article Name	Protein A-coupled Magnetic Beads (recommended for MPCLIA), E. coli
Biozol Catalog Number	ABS-MPC-A001-10MG
Supplier Catalog Number	MPC-A001-10mg
Alternative Catalog Number	ABS-MPC-A001-10MG
Manufacturer	AcroBiosystems
Host	E. coli
Category	Sonstiges
Species Reactivity	Other
Product Description	Protein Acoupled Magnetic Beads (recommended for MPCLIA) is expressed from Recombinant E. coli. - Sonstiges...
Form	Powder
Target	Protein A
Application Notes	The Protein A-coupled Magnetic Beads is used to capture the most antibodies, it can combination with Acridinium ester markers in chemiluminescence technology, The Acridinium ester markers such as Streptavidin-Acridinium ester can capture the biotinylated proteins or molecules, this allows detection of antigen and antibody binding or antibody screening.