

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Anti-FITC Antibody coupled Magnetic Beads (recommended for MPCLIA), Human ABS-MPC-A009-100MG

Article Name	Anti-FITC Antibody coupled Magnetic Beads (recommended for MPCLIA), Human
Biozol Catalog Number	ABS-MPC-A009-100MG
Supplier Catalog Number	MPC-A009-100mg
Alternative Catalog Number	ABS-MPC-A009-100MG
Manufacturer	AcroBiosystems
Host	Human
Category	Sonstiges
Species Reactivity	Rabbit
Product Description	AntiFITC Antibody coupled Magnetic Beads (recommended for MPCLIA) is expressed from Rabbit HEK293. - Sonstiges...
Form	Powder
Target	Anti-FITC
Application Notes	The Anti- FITC Antibody coupled Magnetic Beads is used to capture the FITC Antibody, it can combination with Acridinium markers in chemiluminescence technology, The Acridinium ester markers such as Anti-Human IgG Antibody-Acridinium ester can capture the antibodies , this allows detection of antigen and antibody binding or antibody screening.