

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EasyStep Human VASP(Vasodilator Stimulated Phosphoprotein) ELISA Kit EBT-ELK025ES

Article Name	EasyStep Human VASP(Vasodilator Stimulated Phosphoprotein) ELISA Kit
Biozol Catalog Number	EBT-ELK025ES
Supplier Catalog Number	ELK025ES
Alternative Catalog Number	EBT-ELK025ES-96, EBT-ELK025ES-48, EBT-ELK025ES-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Human
Concentration	10 ng/mL
Range	0.16-10 ng/mL
Sensitivity	17.5 pg/mL
UniProt	P50552
Samples	Serum, plasma

Application Notes	<p>Assay Type: Sandwich. Assay length: 1.5h. Research Area: Signal transduction,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human VASP, and the Human VASP standard plate wells that pre-coated using protein-related techniques are provided separately.</p> <p>Standard/Sample Diluent Buffer or samples are added to the appropriate microtiter plate wells ,then added a HRP-conjugated antibody specific to Human VASP. After TMB substrate solution is added, only those wells that contain Human VASP and HRP-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Human VASP in the samples is then determined by comparing the OD of the samples to the standard curve</p>
-------------------	---