

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Cattle BALP(Bone-specific Alkaline Phosphatase) ELISA Kit EBT-ELK0331

Article Name	Cattle BALP(Bone-specific Alkaline Phosphatase) ELISA Kit
Biozol Catalog Number	EBT-ELK0331
Supplier Catalog Number	ELK0331
Alternative Catalog Number	EBT-ELK0331-96, EBT-ELK0331-48, EBT-ELK0331-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Bovine
Concentration	100 ng/mL
Range	1.57-100 ng/mL
Sensitivity	0.49 ng/mL
Samples	serum, plasma and other biological fluids

Application Notes	<p>Assay Type: Sandwich. Assay length: 3.5h. Research Area: Enzyme &amp; Kinase, Metabolic pathway, Hepatology, Bone metabolism,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle BALP. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle BALP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle BALP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Cattle BALP in the samples is then determined by comparing the OD of the samples to the standard curve</p>
-------------------	---