

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EasyBlot Sampler Kit GTX300111

Article Name	EasyBlot Sampler Kit
Biozol Catalog Number	GTX300111
Supplier Catalog Number	GTX300111
Alternative Catalog Number	GTX300111-1
Manufacturer	GeneTex
Category	Antikörper
Product Description	EasyBlot anti-Rabbit/Mouse IgG is an HRP-conjugated second step reagent that specifically reacts with the native, non-reduced form of rabbit/Mouse IgG and does not bind the reduced, denatured form. Using Easyblot decreases the interference caused by ...
Form	Liquid
Application Notes	Suggested starting dilutions for WB is 1:1000. Use 2-5 µg of the IP antibody, do not exceed 5 µg. Completely reduce all IP samples by 100mM DTT (fresh prepared). For further detail, please refer to the protocol. Not yet tested in other applications. Optimal working dilutions should be determined experimentally by the end user.