

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Precast Gel, Tris-Gly 4-20%, 15 Wells, 1.0mm ABB-BR00001

Article Name	Precast Gel, Tris-Gly 4-20%, 15 Wells, 1.0mm
Biozol Catalog Number	ABB-BR00001
Supplier Catalog Number	BR00001
Alternative Catalog Number	ABB-BR00001-1PCS, ABB-BR00001-10PCS
Manufacturer	ABclonal
Category	Kits/Assays
Product Description	Precast Gel is a safe, fast, and high-performance precast polyacrylamide gel, commonly used in PAGE (Polyacrylamide Gel Electrophoresis) and Western blot detection....
Concentration	Identification of Protein Bands
Application Dilute	NameSpecificationStorage Precast Gel, Tris-Gly 4-20%, 15 Wells, 1.0mm10pcs2-8°C Instruction Manual 1copy /
Application Notes	Cross-Reactivity: 1.Coated plastic gel plates are used to effectively reduce non-specific protein adsorption,resulting in sharper and clearer protein bands. 2.The gel clamps are easy to open,and the plastic plates can be easily pried apart with a screwdriver. 3.Automated gel casting production technology is adopted to ensure high stability and repeatability of product quality. 4.A variety of concentrations of homogeneous gels and gradient gels are available.