

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Tris-Glycine-Native Electrophoresis Buffer Powder ABB-BR00006

Article Name	Tris-Glycine-Native Electrophoresis Buffer Powder
Biozol Catalog Number	ABB-BR00006
Supplier Catalog Number	BR00006
Alternative Catalog Number	ABB-BR00006-10X500ML, ABB-BR00006-5X500ML
Manufacturer	ABclonal
Category	Sonstiges
Product Description	Tris-Glycine-Native Electrophoresis Buffer is compatible with Tris-based gels and is commonly used for separating large-molecular-weight proteins. It is designed for native polyacrylamide gel electrophoresis (Native-PAGE). This product is supplied in ...
Concentration	Protein electrophoresis
Application Dilute	Product CodeProduct NameSpecificationStorage Temperature BR00006Tris-Glycine-Native Electrophoresis Buffer Powder5*500 mL / 10*500 mL Store at room temperature
Application Notes	Cross-Reactivity: 1.Fast. 2.Clear protein bands. 3.Simple operation,no need to weigh or filter reagents.