

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### High Sensitive ELISA Kit for Lipopolysaccharide (LPS) USC-HEB526GE

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| Article Name               | High Sensitive ELISA Kit for Lipopolysaccharide (LPS)                                                |
| Biozol Catalog Number      | USC-HEB526GE                                                                                         |
| Supplier Catalog Number    | HEB526Ge                                                                                             |
| Alternative Catalog Number | USC-HEB526GE-48, USC-HEB526GE-96                                                                     |
| Manufacturer               | Cloud-Clone                                                                                          |
| Category                   | Kits/Assays                                                                                          |
| Species Reactivity         | All                                                                                                  |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                                  |
| Range                      | 0.156-10ng/mL                                                                                        |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 0.058ng/mL                            |
| Samples                    | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Application Notes          | Enzyme-linked immunosorbent assay for Antigen Detection.                                             |