

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Multiplex Assay Kit for Angiotensin II (AngII) ,etc. by FLIA (Flow Luminescence Immunoassay) USC-LMA005RA**

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| Article Name               | Multiplex Assay Kit for Angiotensin II (AngII) ,etc. by FLIA (Flow Luminescence Immunoassay)         |
| Biozol Catalog Number      | USC-LMA005RA                                                                                         |
| Supplier Catalog Number    | LMA005Ra                                                                                             |
| Alternative Catalog Number | USC-LMA005RA-96                                                                                      |
| Manufacturer               | Cloud-Clone                                                                                          |
| Category                   | Kits/Assays                                                                                          |
| Species Reactivity         | Rat                                                                                                  |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                                  |
| Range                      | 1.95-2000pg/mL                                                                                       |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 0.65 pg/mL                            |
| Samples                    | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Application Notes          | FLIA Kit for Antigen Detection.                                                                      |