

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Multiplex Assay Kit for Luteinizing Hormone (LH) ,etc. by FLIA (Flow Luminescence Immunoassay) USC-LMA441RB**

|                            |                                                                                                |
|----------------------------|------------------------------------------------------------------------------------------------|
| Article Name               | Multiplex Assay Kit for Luteinizing Hormone (LH) ,etc. by FLIA (Flow Luminescence Immunoassay) |
| Biozol Catalog Number      | USC-LMA441RB                                                                                   |
| Supplier Catalog Number    | LMA441Rb                                                                                       |
| Alternative Catalog Number | USC-LMA441RB-96                                                                                |
| Manufacturer               | Cloud-Clone                                                                                    |
| Category                   | Kits/Assays                                                                                    |
| Species Reactivity         | Rabbit                                                                                         |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                            |
| Range                      | 29.3-30000pg/mL                                                                                |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 9.767 pg/mL                     |
| Samples                    | Serum, plasma and other biological fluids.                                                     |
| Application Notes          | FLIA Kit for Antigen Detection.                                                                |