

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Multiplex Assay Kit for Somatostatin (SST) ,etc. by FLIA (Flow Luminescence Immunoassay) USC-LMA592MU**

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| Article Name               | Multiplex Assay Kit for Somatostatin (SST) ,etc. by FLIA (Flow Luminescence Immunoassay) |
| Biozol Catalog Number      | USC-LMA592MU                                                                             |
| Supplier Catalog Number    | LMA592Mu                                                                                 |
| Alternative Catalog Number | USC-LMA592MU-96                                                                          |
| Manufacturer               | Cloud-Clone                                                                              |
| Category                   | Kits/Assays                                                                              |
| Species Reactivity         | Mouse                                                                                    |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                      |
| Range                      | 0.49-500pg/mL                                                                            |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 0.163 pg/mL               |
| UniProt                    | <a href="#">P60041</a>                                                                   |
| Samples                    | serum, plasma and other biological fluids                                                |
| Application Notes          | FLIA Kit for Antigen Detection.                                                          |