

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### **Multiplex Assay Kit for Caspase 9 (CASP9) ,etc. by FLIA (Flow Luminescence Immunoassay) USC-LMA627HU**

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| Article Name               | Multiplex Assay Kit for Caspase 9 (CASP9) ,etc. by FLIA (Flow Luminescence Immunoassay) |
| Biozol Catalog Number      | USC-LMA627HU                                                                            |
| Supplier Catalog Number    | LMA627Hu                                                                                |
| Alternative Catalog Number | USC-LMA627HU-96                                                                         |
| Manufacturer               | Cloud-Clone                                                                             |
| Category                   | Kits/Assays                                                                             |
| Species Reactivity         | Human                                                                                   |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                     |
| Range                      | 0.02-20ng/mL                                                                            |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 0.007 ng/mL              |
| UniProt                    | <a href="#">P55211</a>                                                                  |
| Samples                    | Tissue homogenates, cell lysates, cell culture supernates and other biological fluids.  |
| Application Notes          | FLIA Kit for Antigen Detection.                                                         |