

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Multiplex Assay Kit for Kallikrein 7 (KLK7) ,etc. by FLIA (Flow Luminescence Immunoassay) USC-LMB910MU

|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| Article Name               | Multiplex Assay Kit for Kallikrein 7 (KLK7) ,etc. by FLIA (Flow Luminescence Immunoassay) |
| Biozol Catalog Number      | USC-LMB910MU                                                                              |
| Supplier Catalog Number    | LMB910Mu                                                                                  |
| Alternative Catalog Number | USC-LMB910MU-96                                                                           |
| Manufacturer               | Cloud-Clone                                                                               |
| Category                   | Kits/Assays                                                                               |
| Species Reactivity         | Mouse                                                                                     |
| Product Description        | 48T, 96T, 96T*5, 96T*10, 96T*100...                                                       |
| Range                      | 4.88-5000pg/mL                                                                            |
| Sensitivity                | The minimum detectable dose of this kit is typically less than 1.627 pg/mL                |
| UniProt                    | <a href="#">Q91VE3</a>                                                                    |
| Samples                    | Serum, plasma, tissue homogenates and other biological fluids                             |
| Application Notes          | FLIA Kit for Antigen Detection.                                                           |