

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### ExoSparkler EV Protein Labeling Kit - Green DOJ-EX04-10

|                            |                                                                                                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Article Name               | ExoSparkler EV Protein Labeling Kit - Green                                                                                                                                                                                                                   |
| Biozol Catalog Number      | DOJ-EX04-10                                                                                                                                                                                                                                                   |
| Supplier Catalog Number    | EX04-10                                                                                                                                                                                                                                                       |
| Alternative Catalog Number | DOJ-EX04-10                                                                                                                                                                                                                                                   |
| Manufacturer               | Dojindo                                                                                                                                                                                                                                                       |
| Category                   | Kits/Assays                                                                                                                                                                                                                                                   |
| Application                | Exosome Staining                                                                                                                                                                                                                                              |
| Product Description        | Observe exosome dynamics accurately without extracellular aggregation Cover steps from fluorescence labeling to purification<br>*Protein amount : 1-10µg/ sample, Particle count : 10 to 100 x 108 samples (As purified exosome using ultracentrifugation)... |