

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

MDA Assay Kit DOJ-M496-10

| | |
|----------------------------|---|
| Article Name | MDA Assay Kit |
| Biozol Catalog Number | DOJ-M496-10 |
| Supplier Catalog Number | M496-10 |
| Alternative Catalog Number | DOJ-M496-10 |
| Manufacturer | Dojindo |
| Category | Kits/Assays |
| Product Description | Dojindos MDA (malondialdehyde) Assay Kit uses the TBARS method to detect MDA by measuring the fluorescence or absorbance of the adduct of MDA and thiobarbituric acid that develops color according to the amount of MDA. The kit also contains MDA standa... |