

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EasyStep Human NT-ProBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit EBT-ELK035ES

Article Name	EasyStep Human NT-ProBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit
Biozol Catalog Number	EBT-ELK035ES
Supplier Catalog Number	ELK035ES
Alternative Catalog Number	EBT-ELK035ES-96, EBT-ELK035ES-48, EBT-ELK035ES-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Human
Concentration	60 ng/mL
Range	0.94-60 ng/mL
Sensitivity	0.94 ng/mL
Samples	Serum, plasma

Application Notes	Assay Type: Sandwich. Assay length: 1.5h. Research Area: Endocrinology,Cardiovascular biology,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human NT-ProBNP, and the Human NT-ProBNP standard plate wells that pre-coated using protein-related techniques are provided separately. Standard/Sample Diluent Buffer or samples are added to the appropriate microtiter plate wells ,then added a HRP-conjugated antibody specific to Human NT-ProBNP. After TMB substrate solution is added, only those wells that contain Human NT-ProBNP and HRP-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Human NT-ProBNP in the samples is then determined by comparing the OD of the samples to the standard curve
-------------------	--