

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

EasyStep Human SAA(Serum Amyloid A) ELISA Kit EBT-ELK042ES

Article Name	EasyStep Human SAA(Serum Amyloid A) ELISA Kit
Biozol Catalog Number	EBT-ELK042ES
Supplier Catalog Number	ELK042ES
Alternative Catalog Number	EBT-ELK042ES-96, EBT-ELK042ES-48, EBT-ELK042ES-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Human
Concentration	50 ng/mL
Range	0.78-50 ng/mL
Sensitivity	0.15 ng/mL
UniProt	PODJI8
Samples	Serum, plasma

Application Notes

Assay Type: Sandwich. Assay length: 1.5h. Research Area: Infection immunity, Kidney biomarker, Hepatology, Urology,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human SAA, and the Human SAA standard plate wells that pre-coated using protein-related techniques are provided separately. Standard/Sample Diluent Buffer or samples are added to the appropriate microtiter plate wells, then added a HRP-conjugated antibody specific to Human SAA. After TMB substrate solution is added, only those wells that contain Human SAA and HRP-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Human SAA in the samples is then determined by comparing the OD of the samples to the standard curve