

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### EasyStep Human hCG(Hyperglycosylated Chorionic Gonadotropin) ELISA Kit EBT-ELK089ES

Article Name	EasyStep Human hCG(Hyperglycosylated Chorionic Gonadotropin) ELISA Kit
Biozol Catalog Number	EBT-ELK089ES
Supplier Catalog Number	ELK089ES
Alternative Catalog Number	EBT-ELK089ES-96, EBT-ELK089ES-48, EBT-ELK089ES-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Human
Concentration	500mIU/mL
Range	7.81-500 mIU/mL
Sensitivity	0.08 mIU/mL
Samples	Serum, plasma

Application Notes

Assay Type: Sandwich. Assay length: 1.5h. Research Area: Endocrinology, Reproductive science, Neuro science, Genetic science, Hormone metabolism,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Human hCG, and the Human hCG standard plate wells that pre-coated using protein-related techniques are provided separately. Standard/Sample Diluent Buffer or samples are added to the appropriate microtiter plate wells ,then added a HRP-conjugated antibody specific to Human hCG. After TMB substrate solution is added, only those wells that contain Human hCG and HRP-conjugated antibody will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Human hCG in the samples is then determined by comparing the OD of the samples to the standard curve