

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Mouse VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit EBT-ELK10552

Article Name	Mouse VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit
Biozol Catalog Number	EBT-ELK10552
Supplier Catalog Number	ELK10552
Alternative Catalog Number	EBT-ELK10552-96, EBT-ELK10552-48, EBT-ELK10552-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Mouse
Concentration	200 ng/mL
Range	3.13-200 ng/mL
Sensitivity	1.21 ng/mL
UniProt	P55284
Samples	tissue homogenates, cell lysates and other biological fluids

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: CD & Adhesion molecule, Tumor immunity, Infection immunity,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Mouse VE-Cadherin. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Mouse VE-Cadherin. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Mouse VE-Cadherin, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Mouse VE-Cadherin in the samples is then determined by comparing the OD of the samples to the standard curve