

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Cattle SAA(Serum Amyloid A) ELISA Kit EBT-ELK10656

Article Name	Cattle SAA(Serum Amyloid A) ELISA Kit
Biozol Catalog Number	EBT-ELK10656
Supplier Catalog Number	ELK10656
Alternative Catalog Number	EBT-ELK10656-96, EBT-ELK10656-96X5, EBT-ELK10656-48
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Bovine
Concentration	200 ng/mL
Range	3.13-200 ng/mL
Sensitivity	0.273 ng/mL
UniProt	P35541
Samples	serum, plasma and other biological fluids

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Infection immunity, Kidney biomarker, Hepatology, Urology,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle SAA. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle SAA. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle SAA, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Cattle SAA in the samples is then determined by comparing the OD of the samples to the standard curve