

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Dog ELA2(Elastase 2, Neutrophil) ELISA Kit EBT-ELK10700

Article Name	Dog ELA2(Elastase 2, Neutrophil) ELISA Kit
Biozol Catalog Number	EBT-ELK10700
Supplier Catalog Number	ELK10700
Alternative Catalog Number	EBT-ELK10700-96, EBT-ELK10700-96X5, EBT-ELK10700-48
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Canine
Concentration	100 ng/mL
Range	1.57-100 ng/mL
Sensitivity	0.49 ng/mL
Samples	Serum, plasma and other biological fluids

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Enzyme & Kinase, Infection immunity,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Dog ELA2. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Dog ELA2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Dog ELA2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Dog ELA2 in the samples is then determined by comparing the OD of the samples to the standard curve