

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

## Product Datasheet

### Rabbit TIMP4(Tissue Inhibitors Of Metalloproteinase 4) ELISA Kit EBT-ELK2834

Article Name	Rabbit TIMP4(Tissue Inhibitors Of Metalloproteinase 4) ELISA Kit
Biozol Catalog Number	EBT-ELK2834
Supplier Catalog Number	ELK2834
Alternative Catalog Number	EBT-ELK2834-96, EBT-ELK2834-48, EBT-ELK2834-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Rabbit
Concentration	5000 pg/mL
Range	78.13-5000 pg/mL
Sensitivity	26 pg/mL
UniProt	<a href="#">O97591</a>
Samples	Serum, plasma and other biological fluids.

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Tumor immunity, Cardiovascular biology,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Rabbit TIMP4. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rabbit TIMP4. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Rabbit TIMP4, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Rabbit TIMP4 in the samples is then determined by comparing the OD of the samples to the standard curve