

Please note: This document was created automatically and is not a substitute for the manufacturer's original document.

Product Datasheet

Rabbit PPAR-gamma(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit EBT-ELK4611

Article Name	Rabbit PPAR-gamma(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit
Biozol Catalog Number	EBT-ELK4611
Supplier Catalog Number	ELK4611
Alternative Catalog Number	EBT-ELK4611-96, EBT-ELK4611-48, EBT-ELK4611-96X5
Manufacturer	ELK Biotechnology
Category	Kits/Assays
Species Reactivity	Rabbit
Concentration	50 ng/mL
Range	0.79-50 ng/mL
Sensitivity	0.29 ng/mL
UniProt	O19052
Samples	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Application Notes

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Metabolic pathway, Endocrinology, Cardiovascular biology, Developmental science, Bone metabolism, . Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Rabbit PPAR-gamma. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rabbit PPAR-gamma. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Rabbit PPAR-gamma, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Rabbit PPAR-gamma in the samples is then determined by comparing the OD of the samples to the standard curve